



Order Mouth•Mate  
now at our website:

[www.armor-dental.com/  
order-mouth-mate-now](http://www.armor-dental.com/order-mouth-mate-now)

## Eliminates Real Pain and Real Problems for Orthodontic Patients


Maintaining good oral hygiene is important for everyone – and during orthodontic treatment it's even more critical. Why? Braces trap food and plaque easily, which can lead to costly, irreversible problems. For proper healing to progress after a surgical tooth exposure, it is important that the patient doesn't disrupt the area and keep the area clean.

### Keep plaque under control and improve hygiene.

Good oral hygiene is essential, but keeping your teeth clean and gums healthy during orthodontics or exposure surgery can be a major challenge. The main enemy is plaque (biofilm), a mass of bacteria that grows in your mouth that can lead to white spots, tooth decay and gum disease. The food stuck around braces or exposures chain quickly leads to the formation of harmful plaque, which produces harmful acids. Those plaque acids can eventually cause white spots and cavities.

### Protect your new smile with Mouth•Mate!

**Mouth•Mate** is the new at-home preventative product that orthodontists and patients love because it helps keep mouths clean and healthy! With new braces, the inner cheeks can get sore from the wires and brackets, making it difficult and uncomfortable to brush. And patients with exposures can have a challenging time keeping the area clean. **Mouth•Mate** Oval and **Mouth•Mate** Circle are the perfect addition to your hygiene routine!



Protect  
your new  
smile with  
**Mouth•Mate!**

Visit our website: [www.mouth-mate.com](http://www.mouth-mate.com)

**Mouth•Mate** can help make home hygiene easier, reduce plaque, and keep your mouth clean and fresh.



**Mouth•Mate Oval** is the perfect addition to enhance patients' orthodontic hygiene routine. Patients with new braces often experience temporary mouth sores, which they can use **Mouth•Mate** to protect while brushing to promote comfort and healing. **Mouth•Mate Oval** can also retract the cheeks to help the toothbrush clean hard to reach areas where bacteria collects. Better brushing on the gum line reduces plaque build-up, inflammation, and gingivitis, which are two of the most common chronic orthodontic patient conditions.

Product No. AD90002 (Box of 10) - \$49.99



**Mouth•Mate Circle** helps patients with the impacted teeth that are surgically exposed with a bonded bracket and chain. Because the gold chain is tied to braces on an adjacent tooth, patients often have difficulty keeping this area clean. **Mouth•Mate Circle** can help patients avoid both the exposed tooth and the short chain so patients can brush the adjacent teeth next to the exposure while not affecting the surgical area.

Product No. AD90004 (Box of 10) - \$49.99



## Mouth•Mate Benefits

### Patient Benefits

- Establish good hygiene habits early with new braces and exposure surgery.
- Retract cheek to help the toothbrush clean hard to reach areas.
- Prevent cavities and white spots.
- Experience pain-free brushing if common gum sores occur.
- Protect exposure wounds from bonded bracket and chains.

### Orthodontist Benefits

- Improve post hygiene at home for multi-month/year orthodontic patients.
- Reduce plaque, inflammation, and gingivitis.
- Decrease risk of developing white spots and cavities.
- Improve healing and hygiene after a surgical tooth exposure so patient doesn't disrupt the area.
- Eliminate patient's fear associated with brushing with new braces.
- Enhance patient comfort and satisfaction, which improve positive reviews on social media.

**ARMORDENTAL**

We're making it easier for you to help your patients heal better at home.

Armor Dental, whose mission is to modernize patient recovery at home, is a privately owned Boston based manufacturer of dental products that aid in the healing of post surgical wounds and chronic sores. Armor Dental is developing an entirely new category of dental hygiene products to protect oral wounds and canker sores anywhere in the mouth.

**Contact ArmorDental**

20 Camelot Dr.  
Hingham, MA 02043  
206.852.8427

info@armor-dental.com  
www.armor-dental.com